

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST JONATHAN FEARIS AGENCY LA PAZ CO S.O.

DATE 5/22/16 TIME 1438

INTOXILYZER SERIAL # 80-000917 LOCATION QUARTZSITE P.D.

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 099 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: QUARTZSITE PD

Serial Number: 80-000917

Core Version: 8105.48

05/22/2016

14:38:30

Standard Lot#: 0666100A1

Last Changed By: J. FERRIS #608

QAS: J. FERRIS #608

LA HAZ CSO

Operator: J. FERRIS #608

LA HAZ CSO

Test g/210L Time

Air Blank 0.000 14:39:10

0.100 Cal Check 0.999 14:39:31

Air Blank 0.000 14:40:01

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST JONATHAN FRERIS AGENCY LA PAZ CO S.O.
DATE 5/22/16 TIME 1444
INTOXILYZER SERIAL # 80-000917 LOCATION QUARTERSITE P.D.

1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

1. Clock time check.
 2. Date check.

OPERATIONAL TESTS

1. Alcohol-free subject test result 0. 000 AC.
 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
 4. Standard Calibration Check standard 0. 100 AC. Result 0. 197 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

ETHANOL TANK WAS CHANGED PRIOR TO TESTS

SIGNATURE [Signature]

INTOXILYZER 8000
Location: QUARTZSITE PD
Serial Number: 80-000917
Core Version: 8105.48
05/22/2016

14:44:13

Standard Lot#: 06616100A1
Last Changed By: J. FERRIS #608

OAS: J. FERRIS #608
LA PAZ CSO

Operator: J. FERRIS #608
LA PAZ CSO

Subject: SOAP, TEST
DOB: 01/01/1984
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:45:19
Diagnostic Test	Pass	14:45:54
Air Blank	0.000	14:46:22
0.100 Cal Check	0.099	14:46:42
Air Blank	0.000	14:47:12
Subject Test	0.000	14:47:42
Air Blank	0.000	14:48:09
Wait	ABT*	
Air Blank	IBS**	14:48:27

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: QUARTZSITE PD
Serial Number: 80-000917
Core Version: 8105.48
05/22/2016

14:49:12

Standard Lot#: 06616100A1
Last Changed By: J. FERRIS #608

QAS: J. FERRIS #608
LA PAZ CSO

Operator: J. FERRIS #608
LA PAZ CSO

Subject: SOAP, TEST
DOB: 01/01/1984
Sex: M Weight: 200

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	14:50:08
Diagnostic Test	Pass	14:50:44
Air Blank	0.000	14:51:11
0.100 Cal Check	0.100	14:51:32
Air Blank	0.000	14:52:02
Subject Test	0.105*	14:55:09
Air Blank	0.000	14:55:41
Air Blank	0.000	14:56:08
Subject Test	XXX**	14:56:43
Air Blank	RFI***	14:56:50

*Deficient Sample ✓

**Invalid Sample ✓

Begin new deprivation period

***RFI Detect ✓

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST JONATHAN FEARIS AGENCY LA PAZ CO S.O.

DATE 5/22/16 TIME _____

INTOXILYZER SERIAL # 80-000917 LOCATION QUARTERSITE P.D.

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.097 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: QUARTZSITE PD

Serial Number: 80-000917

Core Version: 8105.48

05/22/2016

14:57:41

Standard Lot#: 06616 100A1

Last Changed By: J. FERRIS #608

QAS: J. FERRIS #608

LA PAZ CSQ

Operator: J. FERRIS #608

LA PAZ CSQ

Test	6/210L	Time
Air Blank	0.000	14:58:13
0.100 Cal Check	0.097	14:58:33
Air Blank	0.000	14:59:03

Preventative Maintenance
Performed -- Timer Reset