

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)  
ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST DEPUTY JOHN M. GOMEZ #220 AGENCY LA PAZ COUNTY SHERIFF'S OFFICE

DATE 08 / 22 / 2015 TIME 13:38 Hrs

INTOXILYZER SERIAL # 80 - 001031 LOCATION LPCSO SALOME STATION

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0. 098 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_

INTOXILYZER 8000  
Location: SALOME  
Serial Number: 80-001031  
Core Version: 8105.48  
08/22/2015 13:38:31

-----  
Standard Lot#: 08915100A2  
Last Changed By: J. GOMEZ #220  
-----

QAS: J. GOMEZ #220  
LA PAZ CSO  
-----

Operator: J. GOMEZ #220  
LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	13:39:15
0.100 Cal Check	0.098	13:39:36
Air Blank	0.000	13:40:06

Preventative Maintenance  
Performed -- Timer Reset

