

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST E Schneider AGENCY LPCCSO
DATE 10/16/15 TIME 0116
INTOXILYZER SERIAL # 80-001031 LOCATION Salome Substation

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .104 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.098 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE E Schneider

INTOXILYZER 8000
Location: SALOME
Serial Number: 80-001031
Core Version: 8105.48
10/16/2015 01:16:31

Standard Lot#: 08915100A2
Last Changed By: J. GOMEZ #220

QAS: E. SCHNEIDER #431
LA PAZ CSO

Operator: E. SCHNEIDER #431
LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	01:17:14
0.100 Cal Check	0.098	01:17:34
Air Blank	0.000	01:18:05

Preventative Maintenance
Performed -- Timer Reset