

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)  
ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST DEPUTY JOHN M. GOMEZ #220 AGENCY LA PAZ COUNTY SHERIFF'S OFFICE

DATE 1 / 24 / 2015

TIME 00:36 Hrs

INTOXILYZER SERIAL # 80 - 001196

LOCATION BACK UP / SALOME SUBSTATION

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the

simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

- (✓) 4. Air blank completed.

- (✓) 5. Calibration check completed. Test results 0.097 AC.

- (✓) 6. Air blank completed.

- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_

INTOXILYZER 8000

Location: SALOME SUBSTATION

Serial Number: 80-001196

Core Version: 8105.48

01/24/2015

00:36:02

-----  
Standard Lot#: 08513100A1

Last Changed By: J. GOMEZ #220  
-----

QAS: J. GOMEZ #220

LA PAZ CSO  
-----

Operator: J. GOMEZ #220

LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	00:36:48
0.100 Cal Check	0.097	00:37:09
Air Blank	0.000	00:37:38

Preventative Maintenance  
Performed -- Timer Reset

A large, stylized handwritten signature in black ink, located in the bottom right corner of the page. The signature is cursive and appears to be the name of the person who performed the maintenance.