

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)**  
**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**  
**STANDARD CALIBRATION CHECK PROCEDURE**

**QA SPECIALIST** DEPUTY JOHN M. GOMEZ #220 **AGENCY** LA PAZ COUNTY SHERIFF'S OFFICE

**DATE** 05 / 21 / 2015 **TIME** 10:33 Hrs

**INTOXILYZER SERIAL #** 80 - 001197 **LOCATION** LPCSO BACKUP / LPCSO ADF

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

**OR**

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0. 100 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

**SIGNATURE** \_\_\_\_\_

INTOXILYZER 8000

Location: BACK UP/LPCSO ADF

Serial Number: 80-001197

Core Version: 8105.48

05/21/2015

10:33:39

-----  
Standard Lot#: 14014100A1

Last Changed By: R. ROSADO #557

-----  
OAS: J. GOMEZ #220

LA PAZ CSO

-----  
Operator: J. GOMEZ #220

LA PAZ CSO

-----  
Test                                  g/210L                                  Time  
-----  
Air Blank                                  0.000                                  10:34:21  
0.100 Cal Check                                  0.100                                  10:34:42  
Air Blank                                  0.000                                  10:35:12

Preventative Maintenance  
Performed -- Timer Reset

