

INTOXILYZER 8000

Location: BACK UP/LPCSO ADF

Serial Number: 80-001197

Core Version: 8105.48

06/23/2016

09:48:06

Standard Lot#: 06616100A1

Last Changed By: J. FERRIS #608

QAS: J. FERRIS #608

LA PAZ CSO

Operator: J. FERRIS #608

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	09:48:51
0.100 Cal Check	0.099	09:49:12
Air Blank	0.000	09:49:41

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Dep Ferris AGENCY LPCSO
DATE 6-23-16 TIME 0945
INTOXILYZER SERIAL 80-001197 LOCATION JAIL

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0.099 AC.
6. Air blank completed.
7. Remove printed record. ~~Attach the record~~ to the completed checklist.

SIGNATURE [Signature]