

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST E. Schneider AGENCY LPCSO  
DATE 10/16/15 TIME 0512  
INTOXILYZER SERIAL # 80-001197 LOCATION LPCSO Jail

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10d AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 099 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE E Schneider

INTOXILYZER 8000

Location: BACK UP/LPCSO ADF

Serial Number: 80-001197

Core Version: 8105.48

10/16/2015

05:12:10

-----  
Standard Lot#: 14014100A1

Last Changed By: R. ROSADO #557  
-----

QAS: E. SCHNEIDER #431

LA PAZ CSO  
-----

Operator: E. SCHNEIDER #431

LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	05:12:43
0.100 Cal Check	0.099	05:13:03
Air Blank	0.000	05:13:33

Preventative Maintenance  
Performed -- Timer Reset