

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST E. Schneider AGENCY LPCSO

DATE 4/25/15 TIME 1756

INTOXILYZER SERIAL # 80-001198 LOCATION Ehrenberg Substation

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. .09 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE E Schneider

INTOXILYZER 8000

Location: EHRENBURG STATION

Serial Number: 80-001198

Core Version: 8105.48

04/25/2015

17:56:04

Standard Lot#: 08513100A1

Last Changed By: M. JOHNSON #531

QAS: E. SCHNEIDER #431

LA PAZ CSO

Operator: E. SCHNEIDER #431

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	17:56:37
0.100 Cal Check	0.093	17:56:58
Air Blank	0.000	17:57:27

Preventative Maintenance
Performed -- Timer Reset