

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. ROSADO #557 AGENCY LA PAZ COUNTY SO

DATE 12-24-14 TIME 00:20
INTOXILYZER SERIAL # 80-001198 CAPAZSO LOCATION LPCSO DJ-EMERGENCY

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 094 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE R. ROSADO #557

DPS Form Exh G-2 (Rev 05-01)

Historical Note

New Exhibit G-2 made by final rulemaking at 12 A.A.R. 1916, effective 9:00 a.m., May 18, 2006 (Supp. 06-2).

INTOXILYZER 8000

Location: EHRENBURG STATION

Serial Number: 80-001198

Core Version: 8105.48

12/24/2014

00:20:20

Standard Lot#: 08513100A1

Last Changed By: M. JOHNSON #531

QAS: R. ROSADO #557

LA PAZ CSO

Operator: R. ROSADO #557

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	00:20:49
0.100 Cal Check	0.094	00:21:10
Air Blank	0.000	00:21:39

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4

THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST R. Rosado # 557

AGENCY LA PAZ COUNTY SO

DATE 12-21-14

00:07

TIME

INTOXILYZER SERIAL # 80-001198 LAPAZSO

LOCATION LPCSO D2 - EMBREBERG

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS


- () 1. Clock time check.
- () 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0. 000 AC.
- () 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- () 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- () 4. Standard Calibration Check standard 0. 100 AC. Result 0. 094 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS _____

SIGNATURE 

DPS Form Exh G-4 (Rev 05-01)

Historical Note

New Exhibit G-4 made by final rulemaking at 12 A.A.R. 1916, effective 9:00 a.m., May 18, 2006 (Supp. 06-2).

Operator: R. ROSADO
LA PAZ CSO

#557

Subject: SQAP, TEST,
DOB: 01/01/1990
Sex: M

Weight: 900

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:08:42
Diagnostic Test	Pass	00:09:17
Air Blank	0.000	00:09:45
0.100 Cal Check	0.095	00:10:06
Air Blank	0.000	00:10:35
Subject Test	0.000	00:11:18
Air Blank	ABT*	00:11:24
Air Blank	IPS**	00:11:32

*Sequence Aborted

**Improper Sample

Not a Successfully Completed Test Sequence

Operator: R. ROSADO
LA PAZ CSO

#557

Subject: SQAP,TEST,

DOB: 01/01/1990

Sex: M

Weight: 900

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	00:13:12
Diagnostic Test	Pass	00:13:47
Air Blank	0.000	00:14:14
0.100 Cal Check	0.095	00:14:35
Air Blank	0.000	00:15:04
Subject Test	0.098*	00:18:12
Air Blank	0.000	00:18:44
Air Blank	0.000	00:19:12
Subject Test	XXX**	00:19:55
Air Blank	RFI***	00:20:06

*Deficient Sample

**Invalid Sample

Begin new deprivation period

***RFI Detect

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)


ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. ROSADO # 557 AGENCY LA PAZ COUNTY SO

DATE 10-24-14 TIME 00:05
INTOXILYZER SERIAL # 80-001198 LA PAZ SO LOCATION LPCSO DJ-EURENBERG

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34° C ± 0.2° C
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 095 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE  557

DPS Form Exh G-2 (Rev 05-01)

Historical Note
New Exhibit G-2 made by final rulemaking at 12 A.A.R. 1916, effective 9:00 a.m., May 18, 2006 (Supp. 06-2).

INTOXILYZER 8000

Location: EHRENBURG STATION

Serial Number: 80-001198

Core Version: 8105.48

12/24/2014

00:05:02

Standard Lot#: 08513100A1

Last Changed By: M. JOHNSON #531

QAS: R. ROSADO #557

LA PAZ CSO

Operator: R. ROSADO #557

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	00:05:31
0.100 Cal Check	0.095	00:05:52
Air Blank	0.000	00:06:22

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104(A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST R. ROSADO #557

AGENCY LA PAZ COUNTY SO

DATE 12-24-14

00:05

TIME

INTOXILYZER SERIAL # 80-001198 LAPAZ SO

LOCATION LPCSO DD - FILLCENBERG

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ} \text{C} \pm 0.2^{\circ} \text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 095 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE

R. Rosado #557

DPS Form Exh G-2 (Rev 05-01)

Historical Note

New Exhibit G-2 made by final rulemaking at 12 A.A.R. 1916, effective 9:00 a.m., May 18, 2006 (Supp. 06-2).

INTOXILYZER 8000

Location: EHRENBURG STATION

Serial Number: 80-001198

Core Version: 8105.48

12/24/2014

00:05:02

Standard Lot#: 08513100A1

Last Changed By: M. JOHNSON #531

QAS: R. ROSADO #557

LA PAZ CSO

Operator: R. ROSADO #557

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	00:05:31
0.100 Cal Check	0.095	00:05:52
Air Blank	0.000	00:06:22

Preventative Maintenance
Performed -- Timer Reset