

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Dep Johnson AGENCY LPCSO
DATE 12-16-10 TIME 1610
INTOXILYZER SERIAL # 80-001200 LOCATION COUNTY PARK

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- 72-73-74
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed
5. Calibration check completed. Test results 0. 100 AC.
6. Air blank completed
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000
Location: QTZ 103/LPCS BACK UP
Serial Number: 80-001200
Core Version: 8105.48
12/16/2015 16:10:26

Standard Lot#: 08915100A2
Last Changed By: R. ROSADO #557

QAS: M. JOHNSON #531
LA PAZ CSO

Operator: M. JOHNSON #531
LA PAZ CSO

Test g/210L Time

Air Blank 0.000 16:10:56
0.100 Cal Check 0.100 16:11:17
Air Blank 0.000 16:11:47

Preventative Maintenance
Performed -- Timer Reset