

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST E. Schneider AGENCY LPCS  
DATE 10/16/15 TIME 0351  
INTOXILYZER SERIAL # 80-003861 LOCATION Boat House

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 100 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000  
Location: LPCSO BOATHOUSE  
Serial Number: 80-003861  
Core Version: 8105.48  
10/16/2015 03:51:08

-----  
Standard Lot#: 14014100A1  
Last Changed By: E. SCHNEIDER #431  
-----

QAS: E. SCHNEIDER #431  
LA PAZ CSO  
-----

Operator: E. SCHNEIDER #431  
LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	03:51:48
0.100 Cal Check	0.100	03:52:08
Air Blank	0.000	03:52:38

Preventative Maintenance  
Performed -- Timer Reset